





SUMMER SPECIAL OFFER

1st JULY— 31st AUGUST 2010

Taq DNA Polymerase with Direct Loading Buffer



Taq-LOAD offers a very high PCR sensitivity, leading to greater yields of specific PCR products. Detection of as little as 1fg of DNA has been achieved. Taq DNA Polymerase with Direct Loading buffer offers robust amplification equal, and in some cases superior, to Taq DNA Polymerase with conventional buffers.



**FAST AND
CONVENIENT**

**DIRECT LOADING
ONTO GELS**

**HIGH EFFICIENCY
TAQ DNA POLYMERASE**

EPTAL110 TAQ & LOAD READY TO USE PCR MIX(1000x50µl rxn)

€475

5X concentrated ready to use Mastermix including Taq DNA Polymerase, dNTPs, MgCl₂ and direct loading buffer

Z 100 M CENTRIFUGE

HERMLE
LABORTECHNIK

The small footprint allows the **Z 100 M** to fit on the bench, even there is no space!
Each **Z 100 M** Centrifuge is supplied with 6 x 1,5/ 2,0 ml rotor.
With an on/off switch you can stop the rotor before opening the lid.

Specification:

max. Speed: 6000 rpm
max RCF: 2000 xg
max. Volume: 6 x 1,5/2,0 ml
Dimensions (wxhxd): 15 x 10 x 15



€215

Ref.no.S001 Please quote this reference number on your purchase order to receive this offer.
Prices are subject to VAT and carriage.

The Leinster Technology Centre, Blessington Industrial Estate, Blessington, Co. Wicklow, Ireland
Tel: +353(0)45865440, Fax: +353(0)45865441, Tel (Dublin): +353(0)14962422, Tel (Cork): +353(0)214343551
email: sales@alphatech.ie

The Technology Centre, 2/3 Curran Point, Larne, Co. Antrim, BT40 1AU, N. Ireland
Phone: +44 (0) 28 28260558, Fax: +44 (0) 28 28260548, email: sales@alphatechs.com