





# SUMMER SPECIAL OFFER

1st JULY— 31st AUGUST 2010

## Taq DNA Polymerase with Direct Loading Buffer



Taq-LOAD offers a very high PCR sensitivity, leading to greater yields of specific PCR products. Detection of as little as 1fg of DNA has been achieved. Taq DNA Polymerase with Direct Loading buffer offers robust amplification equal, and in some cases superior, to Taq DNA Polymerase with conventional buffers.



**FAST AND CONVENIENT**

**DIRECT LOADING ONTO GELS**

**HIGH EFFICIENCY TAQ DNA POLYMERASE**

**EPTAL110 TAQ & LOAD READY TO USE PCR MIX(1000x50µl rxn)**

**£370**

5X concentrated ready to use Mastermix including Taq DNA Polymerase, dNTPs, MgCl<sub>2</sub> and direct loading buffer

## Z 100 M CENTRIFUGE

**HERMLE**  
LABORTECHNIK

The small footprint allows the **Z 100 M** to fit on the bench, even there is no space!  
Each **Z 100 M** Centrifuge is supplied with 6 x 1,5/ 2,0 ml rotor.  
With an on/off switch you can stop the rotor before opening the lid.

### Specification:

max. Speed: 6000 rpm  
max RCF: 2000 xg  
max. Volume: 6 x 1,5/2,0 ml  
Dimensions (wxhxd): 15 x 10 x 15



**£174**

Ref.no.S001 Please quote this reference number on your purchase order to receive this offer.  
Prices are subject to VAT and carriage.

The Leinster Technology Centre, Blessington Industrial Estate, Blessington, Co. Wicklow, Ireland  
Tel: +353(0)45865440, Fax: +353(0)45865441, Tel (Dublin): +353(0)14962422, Tel (Cork): +353(0)214343551  
email: sales@alphatech.ie

The Technology Centre, 2/3 Curran Point, Larne, Co. Antrim, BT40 1AU, N. Ireland  
Phone: +44 (0) 28 28260558, Fax: +44 (0) 28 28260548, email: sales@alphatechs.com